

510(k) Summary

ADMINISTRATIVE INFORMATION

| | |
|--------------------------------|---|
| Sponsor/Manufacturer Name | GM dos Reis Indústria e Comércio Ltda Avenida Pierre Simon de La Place 600 Campinas, São Paulo, Brazil 13069-320 Telephone: +55 (19) 3765-9900 Contact person: Paula Oliveira, Quality Manager |
| Submission contact | Graziela Brum Regulatory Affairs Specialist Passarini Regulatory Affairs PR Serviços Regulatórios Administrativos Ltda E-Mail: graziela@rapassarini.com.br Telephone +55 (47) 3804 0075 |
| Preparer | Janine Treter and Graziela Brum Regulatory Affairs Specialist Passarini Regulatory Affairs PR Serviços Regulatórios Administrativos Ltda E-Mail: graziela@rapassarini.com.br Telephone +55 (47) 3804 0075 |
| Date Prepared | 05/march/2021 |
| DEVICE NAME AND CLASSIFICATION | |
| Trade/ Proprietary Name | Versalock Periprosthetic Femur Plates System - GMReis |
| Common Name | Plate, Fixation, Bone; Condylar Plate Fixation Implant; Screw, Fixation, Bone; Cerclage, Fixation. |
| Classification Name | Single/multiple component metallic bone fixation appliances and accessories; Smooth or threaded metallic bone fixation fastener; Bone fixation cerclage. |
| Product Code | HRS, JDP, HWC, JDQ |
| Classification Regulation | 21 CFR 888.3030, 21 CFR 888.3040, 21 CFR 888.3010 |
| Review Panel | Orthopedic |

PREDICATE DEVICE INFORMATION

| | |
|--------------------------|---|
| Primary Predicate Device | K120772 - NCB Periprosthetic Trochanter Plates and Screws and NCB Cable Button for NCB Polyaxial Locking Plate - Zimmer, GmbH |
| Reference Devices | <p>K082527- Zimmer® Universal Locking System: 3.5 mm Locking Plates and Screw (Titanium® Ti-6A1-4V Alloy, CP Grade Titanium) - Zimmer, Inc.</p> <p>K110354- Synthes 4.5mm VA-LCP Curved Condylar Plate System - Synthes (USA)</p> <p>K151907 - Cable-Ready® Cable Grip System: Cable-Ready® GTR System, Cable-Ready® 1.8mm - Zimmer, Inc.</p> <p>K081759 - NCB® Polyaxial Locking Plate System, Proximal Humeral Plates and Zimmer® Universal Locking System; 3.5mm Titanium® Ti-6A1-4V Alloy Locking Screws - Zimmer, GmbH</p> <p>K061211- NCB® Plating System - Zimmer, GmbH</p> <p>K151716- Cable-Ready® Cable Grip System: Cable-Ready® Bone Plate System - Zimmer, Inc.</p> <p>K100111- NCB® Periprosthetic Femur Polyaxial Locking Plate System - Zimmer, Inc.</p> <p>K161616- DePuy Synthes 4.0 mm and 4.5 mm Cortex Screws, DePuy Synthes 2.4 mm Cannulated Screws, DePuy Synthes 3.5 mm and 4.0 mm Cannulated Screws, DePuy Synthes 4.5 mm Cannulated Screws, DePuy Synthes 6.5 mm Cannulated Screws, DePuy Synthes 7.0 mm and 7.3 mm Cannulated Screws, DePuy Synthes 1.5 mm Headless Compression Screws, DePuy Synthes 2.4 mm Headless Compression Screws, DePuy Synthes 3.0 mm Headless Compression Screws, DePuy Synthes 4.5 mm and 6.5 mm Headless Compression Screw - Zimmer, GmbH</p> <p>K042695 Single/multiple component metallic bone fixation appliances and accessories, Smooth threaded metallic bone fixation fastener Zimmer, GmbH</p> <p>K162124 Synthes 4.5mm VA-LCP Curved Condylar Plate System Line Extension, Variable Angle Positioning Pins - Depuy Synthes</p> <p>K182718 - Mini and Micro Fragments Reconstruction System – GMReis - GM dos Reis Industria E Comercio Ltda</p> |

INDICATIONS FOR USE

The Versalock Periprosthetic Femur Plates System – GMReis is indicated for temporary internal fixation and stabilization of fractures and osteotomies of the femur, including:

- Periprosthetic fractures
- Comminuted fractures

- Supracondylar fractures
- Trochanteric fractures
- Fractures in normal and osteopenic bone
- Non-unions and Malunions

DEVICE DESCRIPTION

The purpose of this submission is to obtain marketing clearance for the Versalock Periprosthetic Femur Plates System which is composed of plates, screws and, a cerclage cable and related device.

The subject devices are designed for the treatment of femur fractures, particularly, periprosthetic femur fractures. The plates are available in the following design-types to be used according to the fracture location: Proximal and Distal Femur Periprosthetic Plates, Trochanteric Periprosthetic Plates and Condylar Femur Plates.

The plates are for use with the subject device screws to fix them to the bone. The following compatible screws are available for this purpose: Cortical Screws, Versalock Variable Angle Locking Screws, Versalock Variable Angle Screws, Versalock Variable Angle Cannulated Screws and Versalock Variable Angle Periprosthetic Screw. The Trochanteric Plate Fastening Screw is to connect one plate to another when a Trochanteric Periprosthetic Plates is used. The Versalock Spacer Screw is threaded into the plate hole prior to plate insertion to act as a spacer providing no contact between the plate and the bone surface.

The Gama Cable is a cerclage cable indicated to provide fixation and/or stabilization of the bone when it is not possible the usage of any screw. The Gama Cable related devices are the Gama Cable Lock and the Versalock Connector Screw. During the installation of the Cable, the Gama Cable Lock is crimped to lock the movement of the cable, maintaining the tensioning applied while the Versalock Connector Screw , which holds the cable to the plate and set the proper cable routing position. The Gama Cable and related devices are used in conjunction with the Proximal or Distal Femur Periprosthetic Plates, or Condylar Femur Plates

The subject devices are made of made of titanium alloy (ASTM F136) with exception of the Gama Cable Lock which it is made of commercially pure titanium (ASTM F67). All the subject devices are colored-anodized.

The devices must only be used by qualified surgeons mastering the surgical technique, having been trained and qualified in orthopedic surgeries.

EQUIVALENCE TO MARKETED DEVICE

The subject device is substantially equivalent in indications and design principles to the following predicate device:

K120772 - NCB Periprosthetic Trochanter Plates and Screws and NCB Cable Button for NCB Polyaxial Locking Plate - Zimmer, GmbH

The subject and predicate devices have equivalent intended use and equivalent technological characteristics. The subject and predicate devices are all manufactured from identical or equivalent materials and share equivalent design characteristics. The subject and predicate devices encompass equivalent physical dimensions and are to be sterilized by identical or equivalent method. Any difference in the technological characteristics do not raise new issues of safety or efficacy.

Biocompatibility was established through a risk assessment following ISO 10993-1.

The performance of the subject device plates was demonstrated through static and dynamic testing according to ASTM F382.

The performance of the subject screws was demonstrated through mechanical testing according to ASTM F543.

The performance of the cerclage cable was demonstrated through mechanical testing according to ASTM F2180 and system construct testing.

No clinical data were included in this submission.

The subject devices are provided non-sterile and have no expiration date defined.

CONCLUSION

The documentation submitted in this premarket notification demonstrates that the subject devices have comparable features and performance and, therefore, are substantially equivalent to the identified predicate devices.